

FRIO[®] INSULIN COOLING CASE



No Ice Packs ...

EVER AGAIN!

www.FRIOCase.com

FRÍO[®] INSULIN COOLING CASE

- ▣ *How the FRÍO[®] works*
- ▣ *Sizes and Colors*
- ▣ *Benefits to your customers*

FRÍO[®] INSULIN COOLING CASE

A reusable evaporative cooler...





MEDICATIONS THAT ARE TEMPERATURE SENSITIVE CAN NOT GO OVER 86 F

Most drugs are recommended to be stored at what's known as "controlled room temperature"

FRÍO[®] **INSULIN COOLING CASE**

Exposure to heat could threaten effectiveness of certain medications

INSULIN

VICTOZA

EPI-PENS

MS MEDS

...and many others

BYETTA

BYDUREON

THYROID MEDS

TOPICAL MEDS

An issue for tens of millions of Americans

FRÍO[®] **INSULIN COOLING CASE**

FRIO Cooling Cases

**keep medication at, or below,
room temperature (77-78F)
for up to 2 days
in temperatures of 100F**

FRÍO[®] INSULIN COOLING CASE

Then simply reactivate it
by soaking it again and you are
good for another 2 days,

over and over again...

FRIO[®] INSULIN COOLING CASE



*Sizes and
colors
to meet every
need and
preference...*

FRÍO[®] INSULIN COOLING CASE

- ▣ *Mini – 1 vial (w/syringe)*
- ▣ *Individual– 1 insulin pen*
- ▣ *Duo – 2 vials (w/ syringes) or 2 pens*
- ▣ *Small – 3 vials (w/syringes)*
- ▣ *Large – 6 vials (w/syringes) or 4 pens*
- ▣ *X-Large – 10 vials (w/ syringes) or 8 pens*
- ▣ *Pump – Insulin Pump or 2 vials*

FRIO[®] INSULIN COOLING CASE



Mini

1 vial

FRIO[®] INSULIN COOLING CASE



Individual
1 Pen

FRIO[®] INSULIN COOLING CASE



Duo

2 pens or 2 vials

FRÍO[®] INSULIN COOLING CASE

Small
4 vials



FRÍO[®] INSULIN COOLING CASE

Large
*4 pens or
6 vials*



FRIO[®] INSULIN COOLING CASE



Extra
Large
*8 pens or
10 vials*

FRÍO[®] INSULIN COOLING CASE

Pump

1 pump

or 2 vials



FRÍO[®] INSULIN COOLING CASE

Convenience



FRÍO®



Freedom

FRÍO[®]

INSULIN COOLING CASE

Peace of Mind

FRIO[®] INSULIN COOLING CASE



No Ice Packs ...

EVER AGAIN!

www.FRIOCase.com